

Antibodies Product Datasheet

Product Name:	LRC36 Antibody
Catalog Number:	C44559
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRC36
Concentration:	1 mg/ml
Gene Name:	LRRC36 RORBP70
Protein Name:	Leucine-rich repeat-containing protein 36 (ROR gamma-binding protein 70)
Other Name:	
Human Swiss Prot No.:	Q1X8D7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LRC36
Storage and Stability:	-20°C/1 year