

Antibodies Product Datasheet

Product Name:	LRCH4 Antibody
Catalog Number:	C44568
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRCH4
Concentration:	1 mg/ml
Gene Name:	LRCH4 LRN LRRN1 LRRN4
Protein Name:	Leucine-rich repeat and calponin homology domain-containing protein 4 (Leucine-rich repeat neuronal protein 4) (Leucine-rich neuronal protein)
Other Name:	
Human Swiss Prot No.:	O75427
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LRCH4
Storage and Stability:	-20°C/1 year