

Antibodies Product Datasheet

Product Name:	LRFN1 Antibody
Catalog Number:	C44569
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRFN1
Concentration:	1 mg/ml
Gene Name:	LRFN1 KIAA1484 SALM2
Protein Name:	LRFN1
Other Name:	Leucine-rich repeat and fibronectin type III domain-containing protein 1 (Synaptic adhesion-like molecule 2)
Human Swiss Prot No.:	Q9P244
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.339% sodium azide.
Specificity:	This antibody detects endogenous levels of LRFN1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year