

Antibodies Product Datasheet

Product Name:	LRIF1 Antibody
Catalog Number:	C44571
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRIF1
Concentration:	1 mg/ml
Gene Name:	LRIF1 C1orf103 RIF1
Protein Name:	LRIF1
Other Name:	Ligand-dependent nuclear receptor-interacting factor 1 (Receptor-interacting factor 1)
Human Swiss Prot No.:	Q5T3J3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.328% sodium azide.
Specificity:	This antibody detects endogenous levels of LRIF1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year