

Antibodies Product Datasheet

Product Name:	LRIQ4 Antibody
Catalog Number:	C44575
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRIQ4
Concentration:	1 mg/ml
Gene Name:	LRRIQ4 LRRC64
Protein Name:	LRIQ4
Other Name:	Leucine-rich repeat and IQ domain-containing protein 4 (Leucine-rich repeat-containing protein 64)
Human Swiss Prot No.:	A6NIV6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.75% sodium azide.
Specificity:	This antibody detects endogenous levels of LRIQ4 at Human/Mouse
Storage and Stability:	-20°C/1 year