

Antibodies Product Datasheet

Product Name:	LRP1B Antibody
Catalog Number:	C44578
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 200-280
Concentration:	1 mg/ml
Gene Name:	LRP1B LRPDIT
Protein Name:	Low-density lipoprotein receptor-related protein 1B (LRP-1B) (Low-density lipoprotein receptor-related protein-deleted in tumor) (LRP-DIT)
Other Name:	
Human Swiss Prot No.:	Q9NZR2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LRP1B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year