

Antibodies Product Datasheet

Product Name:	LRRN1 Antibody
Catalog Number:	C44585
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	LRRN1 KIAA1497 Nbla10449 UNQ693/PRO1338
Protein Name:	Leucine-rich repeat neuronal protein 1 (Neuronal leucine-rich repeat protein 1) (NLRR-1)
Other Name:	
Human Swiss Prot No.:	Q6UXK5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LRRN1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year