

Antibodies Product Datasheet

Product Name:	LRSM1 Antibody
Catalog Number:	C44588
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	LRSAM1 TAL UNQ6496/PRO21356
Protein Name:	E3 ubiquitin-protein ligase LRSAM1 (EC 6.3.2.-) (Leucine-rich repeat and sterile alpha motif-containing protein 1) (Tsg101-associated ligase) (hTAL)
Other Name:	
Human Swiss Prot No.:	Q6UWE0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LRSM1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year