

Antibodies Product Datasheet

Product Name:	LS14A Antibody
Catalog Number:	C44589
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 330-410
Concentration:	1 mg/ml
Gene Name:	LSM14A C19orf13 FAM61A RAP55 RAP55A
Protein Name:	Protein LSM14 homolog A (Protein FAM61A) (Protein SCD6 homolog) (Putative alpha-synuclein-binding protein) (AlphaSNBP) (RNA-associated protein 55A) (hRAP55) (hRAP55A)
Other Name:	
Human Swiss Prot No.:	Q8ND56
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LS14A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year