

Antibodies Product Datasheet

Product Name:	LSM2 Antibody
Catalog Number:	C44591
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	LSM2 C6orf28 G7B
Protein Name:	U6 snRNA-associated Sm-like protein LSm2 (Protein G7b) (Small nuclear ribonuclear protein D homolog) (snRNP core Sm-like protein Sm-x5)
Other Name:	
Human Swiss Prot No.:	Q9Y333
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LSM2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year