

Antibodies Product Datasheet

Product Name:	LST8 Antibody
Catalog Number:	C44593
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 100-180
Concentration:	1 mg/ml
Gene Name:	MLST8 GBL LST8
Protein Name:	Target of rapamycin complex subunit LST8 (TORC subunit LST8) (G protein beta subunit-like) (Gable) (Protein GbetaL) (Mammalian lethal with SEC13 protein 8) (mLST8)
Other Name:	
Human Swiss Prot No.:	Q9BVC4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LST8 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year