

Antibodies Product Datasheet

Product Name:	LT4R1 Antibody
Catalog Number:	C44594
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	LTB4R BLT BLT1 BLTR CMKRL1 GPR16 P2RY7
Protein Name:	Leukotriene B4 receptor 1 (LTB4-R 1) (LTB4-R1) (Chemoattractant receptor-like 1) (G-protein coupled receptor 16) (P2Y purinoceptor 7) (P2Y7)
Other Name:	
Human Swiss Prot No.:	Q15722
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LT4R1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year