

Antibodies Product Datasheet

Product Name:	LT4R2 Antibody
Catalog Number:	C44595
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 290-370
Concentration:	1 mg/ml
Gene Name:	LTB4R2 BLT2R BLTR2
Protein Name:	Leukotriene B4 receptor 2 (LTB4-R 2) (LTB4-R2) (LTB4 receptor JULF2) (Leukotriene B4 receptor BLT2) (Seven transmembrane receptor BLTR2)
Other Name:	
Human Swiss Prot No.:	Q9NPC1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	LT4R2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year