

Antibodies Product Datasheet

Product Name:	LTBP2 Antibody
Catalog Number:	C44597
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LTBP2
Concentration:	1 mg/ml
Gene Name:	LTBP2 C14orf141 LTBP3
Protein Name:	Latent-transforming growth factor beta-binding protein 2 (LTBP-2)
Other Name:	
Human Swiss Prot No.:	Q14767
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LTBP2
Storage and Stability:	-20°C/1 year