

Antibodies Product Datasheet

Product Name:	MAP1A Antibody
Catalog Number:	C44671
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1860-1940
Concentration:	1 mg/ml
Gene Name:	MAP1A MAP1L
Protein Name:	Microtubule-associated protein 1A (MAP-1A) (Proliferation-related protein p80) [Cleaved into: MAP1A heavy chain; MAP1 light chain LC2]
Other Name:	
Human Swiss Prot No.:	P78559
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MAP1A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year