

Antibodies Product Datasheet

Product Name:	MARH6 Antibody
Catalog Number:	C44677
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	MARCH6 KIAA0597 RNF176 TEB4
Protein Name:	E3 ubiquitin-protein ligase MARCH6 (EC 6.3.2.-) (Doa10 homolog) (Membrane-associated RING finger protein 6) (Membrane-associated RING-CH protein VI) (MARCH-VI) (Protein TEB-4) (RING finger protein 176)
Other Name:	
Human Swiss Prot No.:	O60337
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MARH6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year