

Antibodies Product Datasheet

Product Name:	MBNL3 Antibody
Catalog Number:	C44692
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MBNL3
Concentration:	1 mg/ml
Gene Name:	MBNL3 CHCR MBLX39 MBXL
Protein Name:	MBNL3
Other Name:	Muscleblind-like protein 3 (Cys3His CCG1-required protein) (Muscleblind-like X-linked protein) (Protein HCHCR)
Human Swiss Prot No.:	Q9NUK0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.58% sodium azide.
Specificity:	This antibody detects endogenous levels of MBNL3 at Human/Mouse
Storage and Stability:	-20°C/1 year