

Antibodies Product Datasheet

Product Name:	MCA3 Antibody
Catalog Number:	C44697
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MCA3
Concentration:	1 mg/ml
Gene Name:	EEF1E1 AIMP3 P18
Protein Name:	Eukaryotic translation elongation factor 1 epsilon-1 (Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3) (Elongation factor p18) (Multisynthase complex auxiliary component p18)
Other Name:	
Human Swiss Prot No.:	O43324
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human MCA3
Storage and Stability:	-20°C/1 year