

Antibodies Product Datasheet

Product Name:	MCCB Antibody
Catalog Number:	C44698
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MCCB
Concentration:	1 mg/ml
Gene Name:	MCCC2 MCCB
Protein Name:	MCCB
Other Name:	Methylcrotonoyl-CoA carboxylase beta chain, mitochondrial (MCCase subunit beta) (EC 6.4.1.4) (3-methylcrotonyl-CoA carboxylase 2) (3-methylcrotonyl-CoA carboxylase non-biotin-containing subunit) (3-methylcrotonyl-CoA:carbon dioxide ligase subunit beta)
Human Swiss Prot No.:	Q9HCC0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.108% sodium azide.
Specificity:	This antibody detects endogenous levels of MCCB at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year