

Antibodies Product Datasheet

Product Name:	ME15P Antibody
Catalog Number:	C44713
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ME15P
Concentration:	1 mg/ml
Gene Name:	METTL15P1 METT5D2
Protein Name:	Putative methyltransferase-like protein 15P1 (EC 2.1.1.-) (Methyltransferase 5 domain-containing protein 2) (Methyltransferase-like protein 15 pseudogene 1)
Other Name:	
Human Swiss Prot No.:	P0C7V9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ME15P
Storage and Stability:	-20°C/1 year