

Antibodies Product Datasheet

Product Name:	MED4 Antibody
Catalog Number:	C44724
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 200-280
Concentration:	1 mg/ml
Gene Name:	MED4 ARC36 DRIP36 VDRIP HSPC126
Protein Name:	Mediator of RNA polymerase II transcription subunit 4 (Activator-recruited cofactor 36 kDa component) (ARC36) (Mediator complex subunit 4) (TRAP/SMCC/PC2 subunit p36 subunit) (Vitamin D3 receptor-interacting protein complex 36 kDa component) (DRIP36)
Other Name:	
Human Swiss Prot No.:	Q9NPJ6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MED4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year