

Antibodies Product Datasheet

Product Name:	MER34 Antibody
Catalog Number:	C44731
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MER34
Concentration:	1 mg/ml
Gene Name:	ERVMER34-1 LP9056
Protein Name:	MER34
Other Name:	HERV-MER_4q12 provirus ancestral Env polyprotein (Endogenous retrovirus group MER34 member 1)
Human Swiss Prot No.:	Q9H9K5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.109% sodium azide.
Specificity:	This antibody detects endogenous levels of MER34 at Human
Storage and Stability:	-20°C/1 year