

Antibodies Product Datasheet

Product Name:	MEX3B Antibody
Catalog Number:	C44740
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MEX3B
Concentration:	1 mg/ml
Gene Name:	MEX3B KIAA2009 RKHD3 RNF195
Protein Name:	MEX3B
Other Name:	RNA-binding protein MEX3B (RING finger and KH domain-containing protein 3) (RING finger protein 195)
Human Swiss Prot No.:	Q6ZN04
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.168% sodium azide.
Specificity:	This antibody detects endogenous levels of MEX3B at Human/Mouse
Storage and Stability:	-20°C/1 year