

Antibodies Product Datasheet

Product Name:	MFSD7 Antibody
Catalog Number:	C44747
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MFSD7
Concentration:	1 mg/ml
Gene Name:	MFSD7 LP2561 UNQ385/PRO717
Protein Name:	MFSD7
Other Name:	Major facilitator superfamily domain-containing protein 7 (Myosin light polypeptide 5 regulatory protein) (MYL5)
Human Swiss Prot No.:	Q6UXD7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.76% sodium azide.
Specificity:	This antibody detects endogenous levels of MFSD7 at Human/Mouse
Storage and Stability:	-20°C/1 year