

Antibodies Product Datasheet

Product Name:	MLL4 Antibody
Catalog Number:	C44773
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1060-1140
Concentration:	1 mg/ml
Gene Name:	WBP7 HRX2 KIAA0304 KMT2B MLL2 MLL4 TRX2
Protein Name:	Histone-lysine N-methyltransferase MLL4 (EC 2.1.1.43) (Lysine N-methyltransferase 2B) (KMT2B) (Myeloid/lymphoid or mixed-lineage leukemia protein 4) (Trithorax homolog 2) (WW domain-binding protein 7) (WBP-7)
Other Name:	
Human Swiss Prot No.:	Q9UMN6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MLL4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year