

Antibodies Product Datasheet

Product Name:	MMRN2 Antibody
Catalog Number:	C44782
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MMRN2
Concentration:	1 mg/ml
Gene Name:	MMRN2 EMILIN3
Protein Name:	MMRN2
Other Name:	Multimerin-2 (EMILIN-3) (Elastin microfibril interface located protein 3) (Elastin microfibril interfacier 3) (EndoGlyx-1 p125/p140 subunit)
Human Swiss Prot No.:	Q9H8L6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.204% sodium azide.
Specificity:	This antibody detects endogenous levels of MMRN2 at Human/Mouse
Storage and Stability:	-20°C/1 year