

Antibodies Product Datasheet

Product Name:	MOGT1 Antibody
Catalog Number:	C44790
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MOGT1
Concentration:	1 mg/ml
Gene Name:	MOGAT1 DC2 DGAT2L1
Protein Name:	2-acylglycerol O-acyltransferase 1 (EC 2.3.1.22) (Acyl-CoA:monoacylglycerol acyltransferase 1) (MGAT1) (Diacylglycerol O-acyltransferase candidate 2) (hDC2) (Diacylglycerol acyltransferase 2-like protein 1) (Monoacylglycerol O-acyltransferase 1)
Other Name:	
Human Swiss Prot No.:	Q96PD6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human MOGT1
Storage and Stability:	-20°C/1 year