

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | MOGT3 Antibody                                                                                                                                                                                                                                                    |
| <b>Catalog Number:</b>        | C44791                                                                                                                                                                                                                                                            |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                                                                                                                                                 |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                                                                                                                                         |
| <b>Dilutions:</b>             | WB 1:500-2000 ELISA 1:5000-20000                                                                                                                                                                                                                                  |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                                                                                                                  |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                                                                                                 |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                             |
| <b>Immunogen:</b>             | Synthesized peptide derived from human protein . at AA range: 10-90                                                                                                                                                                                               |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                                                           |
| <b>Gene Name:</b>             | MOGAT3 DC7 DGAT2L7 UNQ9383/PRO34208                                                                                                                                                                                                                               |
| <b>Protein Name:</b>          | 2-acylglycerol O-acyltransferase 3 (EC 2.3.1.20) (EC 2.3.1.22) (Acyl-CoA:monoacylglycerol acyltransferase 3) (MGAT3) (Diacylglycerol O-acyltransferase candidate 7) (hDC7) (Diacylglycerol acyltransferase 2-like protein 7) (Monoacylglycerol O-acyltransferase) |
| <b>Other Name:</b>            |                                                                                                                                                                                                                                                                   |
| <b>Human Swiss Prot No.:</b>  | Q86VF5                                                                                                                                                                                                                                                            |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                                                                                                                    |
| <b>Specificity:</b>           | MOGT3 Polyclonal Antibody detects endogenous levels of protein.                                                                                                                                                                                                   |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                                                                      |