

Antibodies Product Datasheet

Product Name:	MRCKG Antibody
Catalog Number:	C44812
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1370-1450
Concentration:	1 mg/ml
Gene Name:	CDC42BPG DMPK2
Protein Name:	Serine/threonine-protein kinase MRCK gamma (EC 2.7.11.1) (CDC42-binding protein kinase gamma) (DMPK-like gamma) (Myotonic dystrophy kinase-related CDC42-binding kinase gamma) (MRCK gamma) (MRCKG) (Myotonic dystrophy protein kinase-like gamma) (Myotonic dy
Other Name:	
Human Swiss Prot No.:	Q6DT37
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MRCKG Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year