

Antibodies Product Datasheet

Product Name:	MUC4 Antibody
Catalog Number:	C44851
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1970-2050
Concentration:	1 mg/ml
Gene Name:	MUC4
Protein Name:	Mucin-4 (MUC-4) (Ascites sialoglycoprotein) (ASGP) (Pancreatic adenocarcinoma mucin) (Testis mucin) (Tracheobronchial mucin) [Cleaved into: Mucin-4 alpha chain (Ascites sialoglycoprotein 1) (ASGP-1); Mucin-4 beta chain (Ascites sialoglycoprotein 2) (ASGP-2)]
Other Name:	
Human Swiss Prot No.:	Q99102
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MUC4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year