

Antibodies Product Datasheet

Product Name:	MUC5B Antibody
Catalog Number:	C44852
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1810-1890
Concentration:	1 mg/ml
Gene Name:	MUC5B MUC5
Protein Name:	Mucin-5B (MUC-5B) (Cervical mucin) (High molecular weight salivary mucin MG1) (Mucin-5 subtype B, tracheobronchial) (Sublingual gland mucin)
Other Name:	
Human Swiss Prot No.:	Q9HC84
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MUC5B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year