

Antibodies Product Datasheet

Product Name:	MUL1 Antibody
Catalog Number:	C44854
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	MUL1 C1orf166 GIDE MAPL MULAN RNF218
Protein Name:	Mitochondrial ubiquitin ligase activator of NFKB 1 (EC 6.3.2.-) (E3 SUMO-protein ligase MUL1) (E3 ubiquitin-protein ligase MUL1) (Growth inhibition and death E3 ligase) (Mitochondrial-anchored protein ligase) (MAPL) (Putative NF-kappa-B-activating protein)
Other Name:	
Human Swiss Prot No.:	Q969V5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MUL1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year