

Antibodies Product Datasheet

Product Name:	MXRA5 Antibody
Catalog Number:	C44858
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	MXRA5
Protein Name:	Matrix-remodeling-associated protein 5 (Adhesion protein with leucine-rich repeats and immunoglobulin domains related to perlecan) (Adlican)
Other Name:	
Human Swiss Prot No.:	Q9NR99
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MXRA5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year