

Antibodies Product Datasheet

Product Name:	MYRIP Antibody
Catalog Number:	C44898
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MYRIP
Concentration:	1 mg/ml
Gene Name:	MYRIP SLAC2C
Protein Name:	Rab effector MyRIP (Exophilin-8) (Myosin-VIIa- and Rab-interacting protein) (Synaptotagmin-like protein lacking C2 domains C) (Slac2-c) (Slp homolog lacking C2 domains c)
Other Name:	
Human Swiss Prot No.:	Q8NFW9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human MYRIP
Storage and Stability:	-20°C/1 year