

Antibodies Product Datasheet

Product Name:	NAA25 Antibody
Catalog Number:	C44903
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NAA25
Concentration:	1 mg/ml
Gene Name:	NAA25 C12orf30 MDM20 NAP1
Protein Name:	NAA25
Other Name:	N-alpha-acetyltransferase 25, NatB auxiliary subunit (Mitochondrial distribution and morphology protein 20) (N-terminal acetyltransferase B complex subunit MDM20) (NatB complex subunit MDM20) (N-terminal acetyltransferase B complex subunit NAA25) (p120)
Human Swiss Prot No.:	Q14CX7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.289% sodium azide.
Specificity:	This antibody detects endogenous levels of NAA25 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year