

Antibodies Product Datasheet

Product Name:	NAV3 Antibody
Catalog Number:	C44921
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 430-510
Concentration:	1 mg/ml
Gene Name:	NAV3 KIAA0938 POMFIL1 STEERIN3
Protein Name:	Neuron navigator 3 (Pore membrane and/or filament-interacting-like protein 1) (Steerin-3) (Unc-53 homolog 3) (unc53H3)
Other Name:	
Human Swiss Prot No.:	Q8IVL0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	NAV3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year