

Antibodies Product Datasheet

Product Name:	NDK5 Antibody
Catalog Number:	C44939
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NDK5
Concentration:	1 mg/ml
Gene Name:	NME5
Protein Name:	Nucleoside diphosphate kinase homolog 5 (NDK-H 5) (NDP kinase homolog 5) (Inhibitor of p53-induced apoptosis-beta) (IPIA-beta) (Testis-specific nm23 homolog) (nm23-H5)
Other Name:	
Human Swiss Prot No.:	P56597
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human NDK5
Storage and Stability:	-20°C/1 year