

Antibodies Product Datasheet

Product Name:	NDUAD Antibody
Catalog Number:	C44948
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	NDUFA13 GRIM19 CDA016 CGI-39
Protein Name:	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 13 (Cell death regulatory protein GRIM-19) (Complex I-B16.6) (CI-B16.6) (Gene associated with retinoic and interferon-induced mortality 19 protein) (GRIM-19) (Gene associated with retinoic and IFN
Other Name:	
Human Swiss Prot No.:	Q9P0J0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	NDUAD Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year