

Antibodies Product Datasheet

Product Name:	PABP2 Antibody
Catalog Number:	C45393
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 170-250
Concentration:	1 mg/ml
Gene Name:	PABPN1 PAB2 PABP2
Protein Name:	Polyadenylate-binding protein 2 (PABP-2) (Poly(A)-binding protein 2) (Nuclear poly(A)-binding protein 1) (Poly(A)-binding protein II) (PABII) (Polyadenylate-binding nuclear protein 1)
Other Name:	
Human Swiss Prot No.:	Q86U42
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	PABP2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year