

Antibodies Product Datasheet

Product Name:	PABP4 Antibody
Catalog Number:	C45394
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 190-270
Concentration:	1 mg/ml
Gene Name:	PABPC4 APP1 PABP4
Protein Name:	Polyadenylate-binding protein 4 (PABP-4) (Poly(A)-binding protein 4) (Activated-platelet protein 1) (APP-1) (Inducible poly(A)-binding protein) (iPABP)
Other Name:	
Human Swiss Prot No.:	Q13310
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	PABP4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year