

Antibodies Product Datasheet

Product Name:	14-3-3-pan Antibody
Catalog Number:	C47404
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human 14-3-3-pan
Concentration:	1 mg/ml
Gene Name:	YWHAB
Protein Name:	14-3-3-pan
Other Name:	14-3-3 protein beta/alpha (Protein 1054) (Protein kinase C inhibitor protein 1) (KCIP-1) [Cleaved into: 14-3-3 protein beta/alpha, N-terminally processed]
Human Swiss Prot No.:	P31946/P61981/P27348/P63104/P31947
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat 14-3-3-pan
Storage and Stability:	-20°C/1 year