

Antibodies Product Datasheet

Product Name:	DRA Antibody
Catalog Number:	C47406
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human DRA.
Concentration:	1 mg/ml
Gene Name:	SLC26A3
Protein Name:	Chloride anion exchanger
Other Name:	SLC26A3; DRA; Chloride anion exchanger; Down-regulated in adenoma; Protein DRA; Solute carrier family 26 member 3
Human Swiss Prot No.:	P40879
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DRA Polyclonal Antibody detects endogenous levels of DRA protein.
Storage and Stability:	-20°C/1 year