

Antibodies Product Datasheet

Product Name:	Erbin Antibody
Catalog Number:	C47407
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Erbin Polyclonal Antibody
Concentration:	1 mg/ml
Gene Name:	ERBB2IP ERBIN KIAA1225 LAP2
Protein Name:	Erbin Polyclonal Antibody
Other Name:	Protein LAP2 (Densin-180-like protein) (ErbB2-interacting protein) (Erbin)
Human Swiss Prot No.:	Q96RT1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Erbin Polyclonal Antibody
Storage and Stability:	-20°C/1 year