

Antibodies Product Datasheet

Product Name:	Crystallin- α C Antibody
Catalog Number:	C47415
Reactivity:	Human
Applications:	WB
Dilutions:	WB: 1:500-1000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of Crystallin- α C
Concentration:	
Gene Name:	HSPB8
Protein Name:	Heat shock protein beta-8
Other Name:	HSPB8; CRYAC; E2IG1; HSP22; PP1629; Heat shock protein beta-8; HspB8; Alpha-crystallin C chain; E2-induced gene 1 protein; Protein kinase H11; Small stress protein-like protein HSP22; YM3454
Human Swiss Prot No.:	Q9UJY1
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Crystallin- α C protein.
Storage and Stability:	-20°C/1 year