

Antibodies Product Datasheet

Product Name:	GAPDH Antibody
Catalog Number:	C47418
Reactivity:	Zebrafish
Applications:	WB
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Glyceraldehyde-3-phosphate dehydrogenase
Concentration:	
Gene Name:	GAPDH
Protein Name:	Glyceraldehyde-3-phosphate dehydrogenase
Other Name:	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH; YM3227
Human Swiss Prot No.:	P04406
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects Zebrafish endogenous GAPDH protein.
Storage and Stability:	-20°C/1 year