

Antibodies Product Datasheet

Product Name:	Glucocorticoid receptor Antibody
Catalog Number:	C47430
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-3000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Glucocorticoid receptor
Concentration:	
Gene Name:	NR3C1
Protein Name:	Glucocorticoid receptor
Other Name:	NR3C1; GRL; Glucocorticoid receptor; GR; Nuclear receptor subfamily 3 group C member 1; YM3226
Human Swiss Prot No.:	P04150
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous GR proteins.
Storage and Stability:	-20°C/1 year