

Antibodies Product Datasheet

Product Name:	CIB1 Antibody
Catalog Number:	C47436
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CIB1
Concentration:	1 mg/ml
Gene Name:	CIB1 CIB KIP PRKDCIP
Protein Name:	CIB1
Other Name:	Calcium and integrin-binding protein 1 (CIB) (Calcium- and integrin-binding protein) (CIBP) (Calmyrin) (DNA-PKcs-interacting protein) (Kinase-interacting protein) (KIP) (SNK-interacting protein 2-28) (SIP2-28)
Human Swiss Prot No.:	Q99828
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CIB1
Storage and Stability:	-20°C/1 year