

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                       |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | ATG4c Rabbit Antibody                                                                                                                                                                 |
| <b>Catalog Number:</b>        | C47460                                                                                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                                                     |
| <b>Applications:</b>          | IHC-p, IF(paraffin section)                                                                                                                                                           |
| <b>Dilutions:</b>             | IHC-p 1:50-300                                                                                                                                                                        |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                 |
| <b>Immunogen:</b>             | Recombinant Protein of ATG4c                                                                                                                                                          |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                               |
| <b>Gene Name:</b>             | ATG4C APG4C AUTL1 AUTL3                                                                                                                                                               |
| <b>Protein Name:</b>          | Cysteine protease ATG4C (EC 3.4.22.-) (AUT-like 3 cysteine endopeptidase) (Autophagin-3) (Autophagy-related cysteine endopeptidase 3) (Autophagy-related protein 4 homolog C)         |
| <b>Other Name:</b>            | Cysteine protease ATG4C (EC 3.4.22.-) (AUT-like 3 cysteine endopeptidase) (Autophagin-3) (Autophagy-related cysteine endopeptidase 3) (Autophagy-related protein 4 homolog C); YM3779 |
| <b>Human Swiss Prot No.:</b>  | Q96DT6                                                                                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                               |
| <b>Specificity:</b>           | The antibody detects endogenous ATG4c protein                                                                                                                                         |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                          |