

Antibodies Product Datasheet

Product Name:	C-Jun (Hydroxylated-p244) Antibody
Catalog Number:	C47468
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human C-Jun (Hydroxylated-p244)
Concentration:	1 mg/ml
Gene Name:	JUN
Protein Name:	C-Jun (Hydroxylated-p244)
Other Name:	Transcription factor AP-1 (Activator protein 1) (AP1) (Proto-oncogene c-Jun) (V-jun avian sarcoma virus 17 oncogene homolog) (p39)
Human Swiss Prot No.:	P05412
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat C-Jun (Hydroxylated-p244)
Storage and Stability:	-20°C/1 year