

Antibodies Product Datasheet

Product Name:	HIF-1 β Antibody
Catalog Number:	C47487
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:1000-2000 IHC: 1:200-500
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of HIF-1 β
Concentration:	
Gene Name:	ARNT
Protein Name:	Aryl hydrocarbon receptor nuclear translocator
Other Name:	ARNT; BHLHE2; Aryl hydrocarbon receptor nuclear translocator; ARNT protein; Class E basic helix-loop-helix protein 2; bHLHe2; Dioxin receptor, nuclear translocator; Hypoxia-inducible factor 1-beta; HIF-1-beta; HIF1-beta; YM3370
Human Swiss Prot No.:	P27540
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous HIF-1 β protein.
Storage and Stability:	-20°C/1 year